



Hybridation *in situ* en ME




Cryo0109NOV1071UA Pb.033.tif
NOV1071 UA Pb
Isabelle
Print Mag = 66006x @ 150 mm
500 nm
HV=80kV
Direct Mag = 60000x



Cryo0109NOV1071UA Pb.034.tif
NOV1071 UA Pb
Isabelle
Print Mag = 110010x @ 150 mm
100 nm
HV=80kV
Direct Mag = 100000x



MILLER00N-0509.359.tif
NOV 886 KCI PTA UA
Isabelle
Print Mag = 16501x @ 150 mm
Acquired mai 20, 2009 at 11:37
500 nm
HV=80kV
Direct Mag = 15000x



MILLER00N-0709.122.tif
NOV1071PTA UA
Martin
Print Mag = 44004x @ 150 mm
Acquired Jun 2, 2009 at 4:16 PM
500 nm
HV=80kV
Direct Mag = 40000x